

# Zeiss Upright Switch On/Off Procedures

1. Switch on microscope (wall plug labelled **Microscope**)
2. Switch on camera (wall plug labelled **Camera**)
3. Log in
4. Run Micro-Manager program  
(configure file (\*.cfg) is in C:\temp)

**Need help?** Contact Yan Gu: 07981875417, office: 2c37, ext:6997

1. Close Micro-Manager program
2. Log off
3. Switch off camera (wall plug labelled **Camera**)
4. Switch off microscope (wall plug labelled **Microscope**)
5. Clean the objective (Only when use immersion oil, with lens tissue)
6. Sign off log book